

# Conversion chart

To reproduce the target shade of prosthesis with each product.

	Gumy No.1 	Gumy No.2 	Gumy No.3 	Gumy No.4 
CERAMAGE & CERAMAGE UP	 <ul style="list-style-type: none"> <li><span style="color: lightblue;">■</span> GUM-T</li> <li><span style="color: pink;">■</span> W* : GUM-V* =1:1</li> <li><span style="color: lightcoral;">■</span> GUM-L</li> <li><span style="color: darkred;">■</span> GUM-O:WO=1:4</li> </ul>	 <ul style="list-style-type: none"> <li><span style="color: lightblue;">■</span> GUM-T</li> <li><span style="color: lightcoral;">■</span> GUM-L</li> <li><span style="color: darkred;">■</span> GUM-D</li> <li><span style="color: darkred;">■</span> GUM-O:WO=3:2</li> </ul>	 <ul style="list-style-type: none"> <li><span style="color: lightblue;">■</span> GUM-T</li> <li><span style="color: darkred;">■</span> GUM-Br*</li> <li><span style="color: darkred;">■</span> GUM-D</li> <li><span style="color: darkred;">■</span> GUM-O:WO=4:1</li> </ul>	 <ul style="list-style-type: none"> <li><span style="color: lightblue;">■</span> AM-V</li> <li><span style="color: darkred;">■</span> GUM-D</li> <li><span style="color: darkred;">■</span> GUM-O:GO=4:1</li> </ul>
VINTAGE ZR	Gum-P	Gum-P : Gum-R=2:1	Gum-DP : Gum-R=1:3	Gum-DP
VINTAGE PRO	Gum-LP : Gum-R : Gum-V=2:1:1	Gum-LP : Gum-R=1:2	Gum-DP : Gum-R=2:3	Gum-DP
DIE COLOR WAX	G1 ( Gum Light )	G2 ( Gum Medium )	G3 ( Gum Dark Medium )	G4 ( Gum Dark )

(\* Only available in CERAMAGE UP)